Ambu
Ideas that work for life


Instruction for use


Ambu® **AuraStraight™**


Single Use Laryngeal Mask - Sterile


For use by medical professionals trained in airway management only



Content

Page

English 4
България 5
Česky 6
Dansk 7
Deutsch 8
Ελληνικά 9
Español 10
Eesti 11
Suomi 12
Français 13
Hrvatski 14
Magyar 15
Italiano 16
日本語 17
Lietuviškai 18
Latviski 19
Nederlands 20
Norsk 21
Polski 22
Português 23
Romania 24
Русский 25
Slovenčina 26
Slovenščina 27
Svenska 28
Türkçe 29
中文 30

1. Intended use

The Ambu AuraStraight is intended for use as an alternative to a facemask for achieving and maintaining control of the airway during routine and emergency anaesthetic procedures in fasted patients.

Intended users: "Professionals trained in airway management".

Have a spare Ambu AuraStraight ready and prepared for immediate use.

Contents sterile unless packaging is damaged or opened.

WARNINGS

The user should be familiar with the following warnings prior to use of the AuraStraight.

- Do not use the AuraStraight before precheck: Visual inspection, air leak test of the device. Failure of any test indicates that the device should not be used.
- All Ambu AuraStraight masks must, prior to use, be investigated and inspected for any potential foreign body. The cuff-protector only serves as protecting the Ambu AuraStraight during storage and transportation and must be removed prior to use.
- Do not use the Ambu AuraStraight if there are any bulges on the cuff or if there are any signs of leakage.
- To avoid trauma, do not use excessive force.
- The Ambu AuraStraight is flammable in the presence of lasers and electrocautery equipment.
- For single use only. Use on other patients can cause cross contamination. Do not soak, rinse, or sterilize this device as these procedures may leave harmful residues or cause malfunction of the device. The design and material used are not compatible with conventional cleaning and sterilization procedures.
- The cuff volume or pressure may change in the presence of nitrous oxide, oxygen, or other medical gases.
- The AuraStraight does not protect the trachea or lungs from the risk of aspiration.

CAUTION

- US federal law restricts this device to be sold to or on the order of a physician.
- Ensure that the device is not in any way damaged before use.
- The patency of the AuraStraight should be reconfirmed after any change in the patient's head or neck position.

For further product information see *www.ambu.com*

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Ambu® is registered trademark of Ambu A/S, Denmark.
Ambu A/S is certified according to ISO 13485.
Ambu A/S