beurer
medical

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(GB) Blood glucose test strips GL44, GL50, GL50 evo

Package contents
- 2 containers each with 25 test strips for use with the **Beurer GL44, GL50 and GL50 evo blood glucose monitors**
- This package insert

Note: Before use check that the seals of the containers are intact. If a container seal is damaged: do not use the test strips!

WARNING
- For use in in-vitro diagnostics only (use outside the body).
- Medical carers and others who use this system on several patients must be aware that all products or objects that come into contact with human blood must be handled, even after cleaning, as though they could transfer pathogens.
- Read the package insert and the Instructions for Use for your blood glucose monitor before using the test strips. For reliable results and for complete customer service and user support, and to maintain the manufacturer's warranty, use these test strips only for the specified systems.
- The container with the test strips includes a desiccating agent that could cause skin and eye irritation if inhaled or swallowed. Keep the container away from small children.
- For external use only (IVD). Do not swallow. Test strips can be fatal if swallowed (choking hazard). If a part has been swallowed, seek medical help immediately.
- A test strip must be used only **once** and on **only** one patient.

Proper use
The test strips are intended for measurement of the blood glucose content in human blood (venous or capillary whole blood) with the **Beurer GL44, GL50 and GL50 evo blood glucose monitors** outside the body (IVD).The test strips are suitable for self-testing.

Instructions for Use and Storage
Read carefully the following information on the use and storage of your test strips. All instructions must be followed to ensure that the test strips will give accurate test results.
- Store test strips in a cool, dry place above 2 °C and below 30 °C. Keep test strips away from direct sunlight or heat. Do not store in a refrigerator or freezer.
- Permitted relative air humidity below 90 %.
- Store test strips in their original container – never place them in a different container.
- Close container tightly immediately after removing a test strip.
- Do not use test strips after the use-by date. If expired test strips are used, test results may be inaccurate. The use-by date can be found on the container beside the hourglass icon.
- The test strips can be stored for up to three months after opening the container. Write this use-by date (date of opening + 3 months /) on the label. The storage time is reduced if it extends beyond the use-by date (see date beside the hourglass icon).
- Do not use the test strips if one of the use-by dates (/) has expired.
- The test strips can be touched anywhere with clean, dry hands.
- Use the test strip for measurement immediately on removal from the container.
- Do not bend, cut or change test strips in any way.
- Do not use test strips that have come into contact with liquids for measurements.

How they work
Test strips enable a quantitative measurement of the glucose level in whole blood. When the slot for blood taking comes into contact with a drop of blood, it is automatically filled by simple capillary action. The blood is drawn into by the absorption slot in the test strip and the monitor measures the glucose content of the blood.
The test is based on the measurement of an electric current that is generated by the chemical reaction of the glucose with the enzyme glucose dehydrogenase (Aspergillus oryzae) on the strip. During the reaction, a mediator transports electrons through the electrode surface and so generates a current.
The monitor analyses the current. The current flow is proportional to the glucose content in the blood sample.
The results are shown on the display of the blood glucose monitor. Only a very small volume of blood is required (0.6 microlitre) and the actual measurement takes about five seconds. The test strips can measure glucose values from 20 to 630 mg/dL (1.1 to 35.0 mmol/L).

Checking test strips with the control solution
When you open a new container of test strips, you will need to check that your monitoring system works correctly with them. Use Beurer control solution LEVEL 3 + LEVEL 4 (REF 464.16). If the test result is in the control range, the system is operating correctly and your procedure is correct. The control range for the function test with control solution is printed on the test strip container.
See the Instructions for Use of your monitoring system and the package insert of the control solution for details of testing with the control solution.

Measuring blood glucose
Read the Instructions for Use for your monitor on taking the blood sample, measuring the blood glucose level, and evaluating the test results.

Chemical components of the test-strip sensor:
- FAD-glucose dehydrogenase 6 %
- Potassium ferricyanide 56 %
- Non-reactive components 38 %

Restrictions
- A hematocrit value (proportion of red blood cells) between 30 % and 55 % has no significant influence on the measurement results.
- A very high or very low haematocrit value (proportion of red blood cells) may cause inaccurate measurements. If the haematocrit value is very high (over 55 %), the blood glucose reading may be too low; if the haematocrit value is very low (under 30 %), it may be too high. If you do not know your haematocrit value, ask your doctor.
- Do not use test strips to monitor the blood glucose of newborns.
- Use only fresh capillary whole blood. Do not use serum or plasma.
- Use capillary blood without squeezing the puncture site. If a limb is crushed, the blood is diluted with tissue fluid and this may lead to an incorrect measuring result.
- Do not use the test strips at altitudes above 7010 m.

Disposing of test strips

WARNING
When disposing of test strips, strictly observe the generally applicable precautions for working with blood. Dispose of all blood samples and materials that you or your patients come into contact with correctly to prevent injury and infection of other persons.

Comparison of measured values with laboratory values
Precision
Three lots of the blood glucose test strips have been tested to assess the precision of the GL44 and GL50 evo blood glucose monitors. This includes a repeat assessment using venous blood and a laboratory precision assessment using the control material. The blood glucose content of the venous blood samples ranges from 42.7 to 418.0 mg/dL (2.4 to 23.2 mmol/L) and control material from three concentrations is used.

Results of the repeat precision measurements

Sample	Venous blood		Grand mean value		Pooled standard deviation	Pooled coef-ficient of var-iation (%)
	mg/dL	mmol/L	mg/dL	mmol/L		
1	42.7	2.4	36.0	2.0	2.0	5.6
2	62.0	3.4	59.2	3.3	3.5	5.9
3	120.5	6.7	127.1	7.1	4.1	3.2
4	201.0	11.2	213.8	11.9	6.7	3.1
5	316.5	17.6	329.9	18.3	10.1	3.1
6	418.0	23.2	433.5	24.1	14.5	3.3

Results of the intermediate precision measurement

Sample	Control material		Grand mean value		Pooled standard deviation	Pooled coef-ficient of var-iation (%)
	mg/dL	mmol/L	mg/dL	mmol/L		
1	70.0	3.9	71.3	4.0	1.0	1.4
2	135.6	7.5	136.3	7.6	1.4	1.1
3	351.5	19.5	350.8	19.5	2.8	0.8

System accuracy
Three lots of blood glucose test strips have been tested to assess the system accuracy of the GL44 and GL50 evo blood glucose monitors and to compare it with the reference method in which capillary whole blood concentrations of 34.4 to 442.8 mg/dL (1.9 to 24.6 mmol/L) have been used.

Result of the system accuracy for glucose concentrations <100 mg/dL (<5.55 mmol/L)

Within ±5 mg/dL (±0.28 mmol/L)	Within ±10 mg/dL (±0.56 mmol/L)	Within ±15 mg/dL (±0.83 mmol/L)
55/180 (30.6 %)	111/180 (61.7 %)	175/180 (97.2 %)

Result of the system accuracy for glucose concentrations ≥100 mg/dL (≥5.55 mmol/L)

Within ±5%	Within ±10%	Within ±15%
220/438 (50.2 %)	357/438 (81.5 %)	422/438 (96.3 %)

Results of the system accuracy for combined glucose concentrations between 34.4 mg/dL (1.9 mmol/L) and 442.8 mg/dL (24.6 mmol/L).

Within ±15 mg/dL or ±15% (Within ±0.83 mmol/L or ±15%)
597/618 (96.6 %)

In comparison to the YSI, the GL44, GL50 and GL50 evo met the EN ISO 15197:2015 standard, whereby 95% of the blood glucose values measured must lie within the following ranges: either ±15 mg/dL (±0.83 mmol/L) of the measured average value when using the reference measuring procedure for blood glucose concentrations <100 mg/dL (<5.55 mmol/L) or ±15% for blood glucose concentrations of ≥100 mg/dL (≥5.55 mmol/L). 99% of the individual measured blood glucose values must fall within zones A and B of the Consensus Error Grid (CEG) for diabetes type 1.

Performance evaluation by the user
A study to assess the glucose values of blood samples of capillary blood from the fingertips, which were obtained from 103 individuals that had no special training, produced the following results: 96.7% within ±15 mg/dL (±0.83 mmol/L) and 95.9% within ±15% of the values obtained in the medical laboratory with glucose concentrations of at least 100 mg/dL (5.55 mmol/L).
You will find further details and information regarding blood glucose results and various technologies in generally relevant specialist medical literature.

Influence / Concentration of tested substances	Blood glucose value		50-100 mg/dL (2.8-5.6 mmol/L)	250-350 mg/dL (13.9-19.4 mmol/L)
Acetaminophen	7 mg/dL	(0.46 mmol/L)	6.6 mg/dL (0.37 mmol/L)	4.5%
Ascorbic acid	4 mg/dL	(0.23 mmol/L)	3.3 mg/dL (0.18 mmol/L)	5.1%
Bilirubin	3.3 mg/dL	(0.06 mmol/L)	0.1 mg/dL (0.01 mmol/L)	-1.4%
Cholesterol	400 mg/dL	(10.34 mmol/L)	-6.8 mg/dL (-0.38 mmol/L)	-6.2%
Creatinine	30 mg/dL	(2.65 mmol/L)	0.0 mg/dL (0.00 mmol/L)	-0.1%
Dopamine	2.2 mg/dL	(0.14 mmol/L)	5.0 mg/dL (0.28 mmol/L)	1.0%
EDTA	5.0 mg/dL	(0.17 mmol/L)	-2.0 mg/dL (-0.11 mmol/L)	-2.4%
Galactose	20 mg/dL	(1.11 mmol/L)	-3.1 mg/dL (-0.17 mmol/L)	0.5%
Gentisic acid	7 mg/dL	(0.45 mmol/L)	7.2 mg/dL (0.40 mmol/L)	2.9%
Glutathione	1 mg/dL	(0.03 mmol/L)	-2.6 mg/dL (-0.14 mmol/L)	-3.7%
Haemoglobin	300 mg/dL	(0.05 mmol/L)	-3.1 mg/dL (-0.17 mmol/L)	-2.6%
Heparin	2.1 mg/dL	(0.0018 mmol/L)	-3.0mg/dL (-0.17 mmol/L)	-1.3%
Ibuprofen	50 mg/dL	(2.43 mmol/L)	-2.6 mg/dL (-0.15 mmol/L)	-1.9%
Icodextrin	1094 mg/dL	(0.64–0.78 mmol/L)	-4.17 mg/dL (-0.23 mmol/L)	-2.9%
L-dopa	2 mg/dL	(0.10 mmol/L)	9.3 mg/dL (0.52 mmol/L)	7.9%
Maltose	278 mg/dL	(7.72 mmol/L)	-1.53 mg/dL (-0.09 mmol/L)	-2.6%
Methyldopa	4 mg/dL	(0.19 mmol/L)	7.3 mg/dL (0.41 mmol/L)	0.9%
Pralidoxime iodide	5 mg/dL	(0.14 mmol/L)	1.7 mg/dL (0.09 mmol/L)	-0.1%
Sodium salicylate	40 mg/dL	(2.50 mmol/L)	-3.1 mg/dL (-0.17 mmol/L)	-0.6%
Tolbutamide	100 mg/dL	(3.70 mmol/L)	0.5 mg/dL (0.03 mmol/L)	-0.8%
Tolazamide	2.5 mg/dL	(0.08 mmol/L)	-2.3 mg/dL (-0.13 mmol/L)	1.8%
Triglyceride	800 mg/dL	(9.37 mmol/L)	-7.50 mg/dL (-0.42 mmol/L)	-4.0%
Uric acid	16.5 mg/dL	(0.98 mmol/L)	6.6 mg/dL (0.37 mmol/L)	1.8%
Xylose	9.5 mg/dL	(0.63 mmol/L)	5.6 mg/dL (0.31 mmol/L)	6.6%

The **Beurer GL44, GL50 and GL50 evo blood glucose monitors** are plasma-calibrated.

The test strips are certified in accordance with the following directives and standards:
IVD (98/79/EC) and EN 13640

Customer service address
Please contact our Customer Service department if you have any questions. Please find the service address in the attached address insert.

OUR OBLIGATIONS TO YOU: Our mission is to supply you with high-quality health products and the best customer service. If you are not completely satisfied with this product, please contact our Customer Service department.

Where do you get test strips?
You can get them from your pharmacy without a prescription or wherever the Beurer blood glucose monitor is sold. If you have any more questions on getting test strips please contact Customer Service.

50 test strips REF 464.15

Beurer GmbH, Söflinger Str. 218, 89077 Ulm, Germany
www.beurer-healthguide.com
www.beurer.com

IVD CE 0483

752.254 – 0116_D_GB_F_I_GR Subject to error and change

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